



**Model:**  
Automatic Light Control System

**Description:**  
*Attachment for Automatic Light Control on B-1000 microscopes*

<b>Function</b>	<p><u>Automatic Light Control</u>: the level of light is adjusted by the microscope in order to maintain the same level as the one the user has chosen, no matter if the aperture of the diaphragm changes, another objective is inserted, opacity of the sample changes, etc.</p> <p>Controlled by the main keyboard on B-1000 body.</p>
<b>Dimensions</b>	<p>Height: 52.5 mm Diameter: 100 mm</p>
<b>NOTES</b>	<p>M-1030 is compatible only with B-1000BF and B-1000PH models.</p>

Όραση οράσκει

 **New York  
Microscope  
Company**  
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
Tel: (877) 877-7274 | Fax: (516) 801-2046  
Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)